

## PsL CLEA

- ***Pseudomonas stutzeri* lipase**
- **Immobilized as Cross-Linked Enzyme Aggregate**



***Pseudomonas stutzeri* lipase (PsL) is a lipase widely available by over expression in a production organism. Lipases catalyze the hydrolysis of triglycerides. The proprietary CLEA methodology has been applied to immobilize this enzyme.**

### **CLEA Methodology**

**Our proprietary methodology to immobilize enzymes as Cross-Linked Enzyme Aggregates (CLEAs) consists of covalent cross linking of precipitated enzymes. This efficient and economically attractive method yields immobilized biocatalysts that do not include support material and therefore have a very high activity per unit volume.**

## Product Properties

<b>Product Type:</b>	Immobilized form of lipase from <i>Pseudomonas stutzeri</i> as a Cross-Linked Enzyme Aggregate (CLEA)
<b>Formulation:</b>	Suspension buffer
<b>Enzyme Type:</b>	Lipase, Triacylglycerol hydrolase, EC 3.1.1.3
<b>Natural Reaction:</b>	Hydrolysis of fats and oils
<b>Substrate Specificity:</b>	Preferentially hydrolysis 1,3 position
<b>Typical activity:</b>	CLEA-ST: 5000 Units/ml <sup>*</sup> ; CLEA-OM: 3500 units/ml <sup>*</sup>

\* 1 unit will catalyse the formation of 1µmol butyric acid from tributyrin at 40°C and pH 7.5

CLEA® is a registered trademark of CLEA Technologies BV.

# *Specific Product Specification*

## *PsL CLEA*

### *Applications*

Lipases in general are used in a wide variety of applications in the fine chemistry, laundry and food industry. In organic synthesis they are used in the production of enantiopure alcohols, amines or acids via ester hydrolysis in aqueous media or via direct esterification in organic media.

### *Storage and Stability*

The PsL CLEA® is best stored in a cool and dry environment. Storage at 4 °C is recommended. Under these conditions the PsL CLEA® retains its activity for at least 12 months. It is recommended to keep the CLEA in suspension upon storage.

### *Formulations*

PsL CLEA is available in two formulations. The standard CLEA and a formulation that is optimized for use in organic media.

**PSL CLEA-ST:**  
Standard formulation of the CLEA

**PSL CLEA-OM:**  
Formulation optimized for use in organic media

### *Pricing and Availability*

PsL CLEA is available in two formulations that are described elsewhere. The available quantities range from 10 kU to giga unit scale. Please inquire for availability, lead times and prices.

### *References*

1. Sheldon, Roger A; Sorgedragger, Menno; Janssen, Michiel H. A. **Use of Cross-linked Enzyme aggregates (CLEAs) for performing biotransformations.** *Chimica oggi, Chemistry Today* 2007, 25(1), 48-52.
2. Sheldon, R. A; Schoevaart, R; Van Langen, L.M. **Cross-linked enzyme aggregates (CLEAs): A novel and versatile method for enzyme immobilization (a review).** *Biocatalysis and Biotransformation* 2005, 23(3/4), 141-147.
3. Sheldon, Roger A; Schoevaart, R; van Langen, Luuk M. **CLEAs: An effective technique for enzyme immobilization.** *Specialty Chem.* 2003, July/August, 40-42.
4. Cao, Linqiu; van Langen, Luuk; Sheldon, Roger A. **Immobilised enzymes: carrier-bound or carrier-free?** *Curr. Opin. Biotechnol.* 2003, 14, 387-394.