

BS protease CLEA

- **Bacillus subtilis protease**
- **Optimum pH 7; optimum temperature 55 °C**
- **Immobilized as Cross-Linked Enzyme Aggregate**



Proteases catalyse the hydrolysis of peptide bonds or are used in peptide synthesis. The proprietary CLEA[®] methodology has been applied to create a highly active immobilised preparation of this protease that is stable against auto-proteolysis, common for proteases.

CLEA Methodology

Our proprietary methodology to immobilize enzymes as Cross-Linked Enzyme Aggregates (CLEAs) consists of covalent cross linking of precipitated enzymes. This efficient and economically attractive method yields immobilized biocatalysts that do not include support material and therefore have a very high activity per unit volume.

Product Properties

Product Type:	Immobilized form of a protease from <i>Bacillus subtilis</i> . Immobilized as a Cross-Linked Enzyme Aggregate (CLEA).
Formulation:	Suspension in buffer
Enzyme Type:	Subtilisin, serine endoprotease EC 3.4.21.62
Natural Reaction:	Peptide synthesis and breakdown
Substrate Specificity:	Non-specific. Neutral or acidic amino acids preferred
Typical activity:	CLEA-ST: ~ 5 Units/ml*

* 1 unit will catalyse the formation of 1µmol N-acetylglycine acid from N-acetylglycine ethyl ester at 40°C and pH 7.5

CLEA[®] is a registered trademark of CLEA Technologies BV.

Specific Product Specification

B. subtilis protease CLEA

Applications

Proteases are used in a wide variety of applications in the laundry and food industry. In organic synthesis they are used in peptide synthesis and for the production of enantiopure alcohols, amines or acids via ester synthesis or hydrolysis.

Storage and Stability

The BS protease CLEA® is best stored in a cool and dry environment. Storage at 4 °C is recommended. Under these conditions the CLEA® retains its activity for at least 12 months.

Formulations

BS CLEA is available as a suspension in buffer in three formulations. The standard CLEA and a formulation that is optimized for use in organic media. Besides a formulation with ultra fine particles is available as CLEA-UF, which has shown itself very efficient in surface treatment of for instance wool fibers.

CLEA-ST:

Standard formulation of Alcalase CLEA

CLEA-OM:

Formulation optimized for use in organic media

CLEA-UF:

CLEAs with a ultra fine particle size

Pricing and Availability

BS protease CLEA is available in the formulations that are described before. The available quantities range from 10 kU to giga unit scale. Please inquire for availability, lead times and prices.

References

1. Sheldon, Roger A; Sorgedragger, Menno; Janssen, Michiel H. A. **Use of Cross-linked Enzyme aggregates (CLEAs) for performing biotransformations.** Chimica oggi, Chemistry Today 2007, 25(1), 48-52.
2. Sheldon, R. A; Schoevaart, R; Van Langen, L.M. **Cross-linked enzyme aggregates (CLEAs): A novel and versatile method for enzyme immobilization (a review).** Biocatalysis and Biotransformation 2005, 23(3/4), 141-147.
3. Sheldon, Roger A; Schoevaart, R; van Langen, Luuk M. **CLEAs: An effective technique for enzyme immobilization.** Specialty Chem. 2003, July/August, 40-42.
4. Cao, Linqiu; van Langen, Luuk; Sheldon, Roger A. **Immobilised enzymes: carrier-bound or carrier-free?** Curr. Opin. Biotechnol. 2003, 14, 387-394.